TEG[®] 6s

Hemostasis System

Citrated: K, KH, RT, FF

REF 07-601

HAEMONETICS[®]
THE Blood Management Company[®]

CE

IVD

Rx Only

HAEMONETICS S.A.
Signy Centre
Rue des Fléchères 6
1274 Signy-Centre
Switzerland

© 2015, 2016 Haemonetics Corporation.

Haemonetics, *THE* Blood Management Company, TEG, RapidTEG, and Thrombelastograph are trademarks or registered trademarks of Haemonetics Corporation in the United States and/or other countries.

PlateletMapping is a registered trademark of Cora Healthcare, Inc.
US Patent Nos. 7,261,861, 7,879,615, 8,236,568, 6,787,363, 9,068,966

EN - English

Proprietary Name and Established Name

Proprietary Name: TEG® 6s System and the following citrated assays and reagents:

- CK (Citrated Kaolin) (Kaolin + CaCl₂, for citrated blood) reagent
- CRT (Citrated RapidTEG) (RapidTEG® + CaCl₂, for citrated blood) reagent
- CKH (Citrated Kaolin with Heparinase) (Kaolin + CaCl₂ + Heparinase, for citrated blood) reagent
- CFF (Citrated Functional Fibrinogen) (Tissue Factor (TF) + CaCl₂ + ReoPro®, for citrated blood) reagent

Established Name: Whole Blood Hemostasis System

Intended Use

The TEG® 6s Hemostasis Analyzer is a non-invasive diagnostic instrument designed to monitor and analyze the coagulation state of a blood sample in order to assist in the assessment of patient clinical hemostasis conditions. The TEG 6s analyzer is indicated for use with adult patients where an evaluation of blood coagulation properties is desired.

Results from the TEG 6s analyzer should not be the sole basis for a patient diagnosis, but should be evaluated together with the patient's medical history, the clinical picture and, if necessary, other coagulation tests.

Summary and Explanation of the Test

The TEG System and the Citrated Multichannel Cartridge are used to test the hemostasis properties of citrated blood samples using four different assays/reagents simultaneously, one in each of the four cartridge channels. For the CK assay, Kaolin-activated test methods are used to reduce variability and to reduce the running time of a native whole blood sample. The CRT assay maximally accelerates the clotting process by simultaneously activating the intrinsic and extrinsic coagulation pathways, allowing for the maximum clot strength to be reached more quickly. CKH is similar to CK but with the addition of heparinase, to neutralize the effects of heparin in the blood. If there is heparin present, the R time for CK will be significantly longer than the R time for CKH. The CFF assay fully inhibits all the platelets, thereby excluding their contribution to clot strength, and therefore measures only the functional fibrinogen contribution to clot strength.

Principles of the Procedure

CK Principle

Kaolin is a standardized reagent consisting of dry Kaolin and 0.85% Saline solution. Kaolin-activated test methods are used to reduce variability and to reduce the running time of a native whole blood sample. Use of these particles of hydrated aluminum silicate shortens coagulation time because Kaolin acts as a contact surface activator (intrinsic pathway), which activates Factor XII and platelets and stimulates the reserve clotting ability of a blood sample. Kaolin is combined with CaCl₂ to neutralize any sodium citrate in the blood.

CRT Principle

RapidTEG maximally accelerates the clotting process by simultaneously activating the intrinsic and extrinsic coagulation pathways using a high concentration of Kaolin and Tissue Factor (TF). This closely reflects the physiological clotting process and yields results significantly faster than a native, Kaolin, and/or TF activated test. CaCl_2 is included to neutralize any sodium citrate in the blood.

CKH Principle

Heparin is commonly used as an anticoagulant in surgical procedures. Even very low concentrations of Heparin, fractions of IU/mL of blood, can noticeably increase the R time and can even completely anticoagulate the blood, making it difficult if not impossible to monitor developing coagulopathies that are masked by high levels of therapeutic Heparin. The Kaolin with Heparinase assay will rapidly and specifically neutralize the anticoagulant property of Heparin. CaCl_2 is included to neutralize any sodium citrate in the blood.

CFF Principle

The Functional Fibrinogen reagent activates the extrinsic pathway using tissue factor and inhibits platelet aggregation using a platelet inhibitor that binds to GPIIb/IIIa receptors. The Functional Fibrinogen reagent fully inhibits all the platelets, thereby excluding their contribution to clot strength (MA), and therefore measures only the functional fibrinogen contribution to clot strength. Functional Fibrinogen is combined with CaCl_2 to neutralize any sodium citrate in the blood. Determining the independent contribution of platelets and fibrinogen further improves the use of the TEG System in helping to detect and diagnose hemostasis defects.

Reagents

CK (Citratd Kaolin) (Kaolin + CaCl_2 , for citrated blood) reagent:

Dried Kaolin and Calcium Chloride (adequate to reverse citration of the sample).

CRT (Citratd RapidTEG) (RapidTEG[®] + CaCl_2 , for citrated blood) reagent:

Dried Kaolin, Tissue Factor and Calcium Chloride (adequate to reverse citration of the sample).

CKH (Citratd Kaolin with Heparinase) (Kaolin + CaCl_2 + Heparinase, for citrated blood) reagent:

Dried Kaolin, Calcium Chloride (adequate to reverse citration of the sample) and Heparinase (adequate to reverse 6 IU heparin/mL of blood).

CFF (Citratd Functional Fibrinogen) (Tissue Factor (TF) + CaCl_2 + ReoPro[®], for citrated blood) reagent:

Dried Tissue Factor, Calcium Chloride (adequate to reverse citration of the sample) and ReoPro[®].

Warnings and Precautions

For In-Vitro Diagnostic Use

Sources of Reagent Error

- Cartridges should be refrigerated (2° C to 8° C) in sealed foil pouches when not in use.
- Patient samples should be tested following a minimum of 15 minutes of incubation after collection in citrated tubes. Samples should not be tested following more than 2 hours post-collection.

Sources of Procedural Error

- Use only 3.2% sodium citrate tubes, completely filled by vacuum. Do not collect blood samples into heparin tubes or any other tubes not specifically recommended.

Sample Preparation

No reconstitution, mixing or dilution of the sample is required.

Storage instructions

Citrated Multichannel cartridges must be protected from temperature extremes above 35°C and must be stored refrigerated at 2°C to 8°C as soon as possible after receipt and prior to use. Cartridges are for single use only and must be discarded after use.

Analyzer

Refer to the *TEG® 6s User Manual*.

Specimen Collection and Preparation for Analysis

For blood draws, use only 3.2% sodium citrate tubes, completely filled by vacuum. No additives or preservatives are necessary to maintain the integrity of the sample, but samples must be used after a 15-minute incubation and within 2 hours of draw to maintain in-use stability.

For sample collection, follow standard phlebotomy procedures (reference CLSI H58-A). After collection, gently invert the tube 5 times and let rest for 15 minutes at room temperature prior to use.

Known Interfering Substances

CK assay was tested for interference factors according to CLSI EP7-A2. Potential interferents tested were Absence of a Discard Tube, Short Draw, Hemolysis, and epsilon aminocaproic acid (EACA). Only Hemolysis was found to be an interferent.

CRT assay was tested for interference factors according to CLSI EP7-A2. Potential interferents tested were Absence of a Discard Tube, Short Draw, Hemolysis, and epsilon aminocaproic acid (EACA). Short Draw and Hemolysis were found to be interferents.

CKH assay was tested for interference factors according to CLSI EP7-A2. The potential interferent tested was Protamine. Results show that Protamine is an interferent at concentrations above 0.062 mg/mL.

CFF assay was tested for interference factors according to CLSI EP7-A2. The potential interferent tested was Heparin. Results show that Heparin is an interferent for Functional Fibrinogen above Heparin concentrations of 1 IU/mL.

Materials Provided

- Ten (10) Citrated Multichannel cartridges, each sealed with a desiccant pack in a foil pouch. Each cartridge contains the following dried reagents: Kaolin, RapidTEG, Kaolin with Heparinase, and Functional Fibrinogen, all with calcium chloride.
- Twelve (12) transfer pipettes

Materials Required (but not supplied)

- TEG® 6s Analyzer (Source: Haemonetics) – see User Manual for instructions for use.
- Blood drawing supplies (any medical supplier)

Test Procedure

1. Remove one cartridge pouch from refrigerated storage.
2. Tear open the pouch, remove the cartridge and discard the pouch and desiccant pack.
3. From the *Home* screen on the TEG analyzer, select **new test**.
4. Select the patient ID you wish to use. If the patient is not listed, you may add a new patient by selecting the + option and entering a patient ID on the screen. Then press **next**.
5. Insert the cartridge into the slot, as indicated on the screen, with the bar code on the left side.
6. After the cartridge pretest has completed and you have verified that the assay is what you intended to run, press **next**. If the cartridge is not the desired assay, press **stop** and start over.
7. On the *Test Information* screen, enter information for the test (if desired), and then press **next**.
8. Using a pipette, load blood into the cartridge sample port, filling up to or above the line indicated on the cartridge.
9. Press **next**.
The TEG analyzer starts the test, and the results are displayed as they become available.
10. When the analyzer displays the "Remove cartridge" prompt, remove the used cartridge from the slot and immediately dispose of it.

Calibration

The TEG analyzer does not require routine calibration.

Quality Control (QC) Procedures

The TEG analyzer performs internal QC checks during a pretest when the cartridge is inserted as well as throughout the duration of the test. This verifies that all electromechanical and pneumatic functions of the analyzer-cartridge combination are operating satisfactorily. Laboratories should perform QC testing as required in accordance with current CLIA regulations regarding levels and frequency of control materials. Haemonetics provides Quality Control (QC) Level 1 and Level 2 materials and Calcium Chloride cartridges to be used for this purpose.

Normal blood and Haemonetics Abnormal QC material may be used with the Citrated Multichannel cartridge to verify the performance of the citrated reagents. This may be performed by the laboratory to verify appropriate handling and storage conditions upon receipt of new cartridges, or whenever cartridges are exposed to adverse environmental conditions. Refer to the Abnormal QC Vial Test Kit (07-662) package insert for details on performing this test.

Results

No calculations are necessary to obtain test results; these are displayed on the TEG analyzer display, and may be printed if a printer is attached.

CK Test Results

The R, K, Angle, MA and LY30 parameters are measured and displayed at the conclusion of the test, along with the corresponding reference ranges established for these parameters.

CRT Test Results

The TEG-ACT, R, K, Angle, MA and LY30 parameters are measured and displayed at the conclusion of the test, along with the corresponding reference ranges established for these parameters.

CKH Test Results

The R, K, Angle, and MA parameters are measured and displayed at the conclusion of the test, along with the corresponding reference range established for these parameters.

CFF Test Results

The MA parameter and the estimated Functional Fibrinogen Level (FLEV) are measured and displayed at the conclusion of the test. FLEV is calculated through a transformation of the MA value.

Limitations of the Procedure

Use of Results

TEG analyzer results should always be considered within the clinical context of the individual patient's case. In the event of inconsistencies with the patient's clinical status, samples should be repeated or supplemented with additional clinical information and/or data including performance verification in the user's patient population.

Expected Values

Expected values for test results are within the Reference Ranges for a reference population that were established according to CLSI C28-A3c. The following tables contain the reference range data for each reagent and parameter. It is recommended that customers validate reference ranges based on their local population.

CK Reference Ranges

Citrated Blood Parameter	N	Range
R (min)	157	4.6 – 9.1
K (min)	157	0.8 – 2.1
Angle (deg)	155	63 – 78
MA (mm)	151	52 – 69
LY30 (%)	132	0.0 – 2.6

CRT Reference Ranges

Citrated Blood Parameter	N	Range
TEG-ACT (sec)	157	82 – 152
R (min)	157	0.3 – 1.1
K (min)	156	0.8 – 2.7
Angle (deg)	154	60 – 78
MA (mm)	152	52 – 70
LY30 (%)	131	0.0 – 2.2

CKH Reference Ranges

Citrated Blood Parameter	N	Range
R (min)	155	4.3 – 8.3
K (min)	157	0.8 – 1.9
Angle (deg)	154	64 – 77
MA (mm)	154	52 – 69

CFF Reference Ranges

Citrated Blood Parameter	N	Range
MA (mm)	151	15 – 32
FLEV (mg/dl)	152	278 – 581

Specific Performance Characteristics

CK Precision

The precision of the Kaolin test was evaluated according to CLSI EP5-A2. Testing was performed for within-run (using three donors), within-device (20-day tests using QC Levels 1 and 2), between lots (using three different cartridge lots), and between-operator (using five operators). Also, multi-site reproducibility tests were performed across three sites, three lots, three operators and three analyzers. Results included Coefficient of Variance (CV) values for all precision tests. The CV of the test results for the Kaolin Assay was less than 10% for all parameters.

CK Sensitivity and Specificity


The Kaolin reagent is often used in conjunction with tests performed on blood samples from patients administered with heparin, with the goal of either determining heparin effect (titration of dose) or determining full reversal (post-protamine) of heparin administered. It is, therefore, important that the Clinical Sensitivity and Specificity for the presence of heparin be determined. Five different heparin spiked samples were examined, with concentrations spread across the heparin therapeutic ranges (0.2 to 6.0 IU/mL). 0.2 IU/mL is the lowest therapeutic dose typically used and 6.0 IU/mL represents a typical dose of heparin used when patients are fully heparinized during cardio-pulmonary by-pass. In addition, five (5) non-heparin spiked samples were analyzed. Kaolin's sensitivity to Heparin was 100%, with the R parameter elongated and out of range for all concentrations. Specificity was also 100%, with all un-spiked samples having an R parameter within normal range.

CRT Precision

The precision of the RapidTEG test was evaluated according to CLSI EP5-A2. Testing was performed for within-run (using three donors), within-device (20-day tests using QC Levels 1 and 2), between lots (using three different cartridge lots), and between-operator (using five operators). Also, multi-site reproducibility tests were performed across three sites, three lots, three operators and three analyzers. Results included Coefficient of Variance (CV) values for all precision tests. The CV of the test results for the RapidTEG Assay was less than 13% for the TEG-ACT parameter and less than 10% for the K, Angle and MA parameters.

CRT Linearity

RapidTEG was tested for linearity according to CLSI EP6-A. The TEG-ACT parameter was found to be linear with respect to concentrations of Heparin, as shown below.

CKH Precision

The precision of the Kaolin with Heparinase test was evaluated according to CLSI EP5-A2. Testing was performed for within-run (using three donors), within-device (20-day tests using QC Levels 1 and 2), between lots (using three different cartridge lots), and between-operator (using five operators). Also, multi-site reproducibility tests were performed across three sites, three lots, three operators and three analyzers. Results included Coefficient of Variance (CV) values for all precision tests. The CV for the test results for the Kaolin with Heparinase Assay was less than 10% for the R parameter.

CKH Neutralization


The Kaolin with Heparinase assay was tested to verify that 6 IU of Heparinase per 1 mL of blood will effectively neutralize the effects of both Unfractionated Heparin (UFH) and Low Molecular Weight Heparin (LMWH) on Kaolin-activated citrated whole blood samples. This is the amount of Heparinase used in the Kaolin with Heparinase assay. Tests were performed using UFH and LMWH spiked samples and un-spiked samples that confirmed the neutralization effectiveness of the assay.

CFF Precision

The precision of the Functional Fibrinogen test was evaluated according to CLSI EP5-A2. Testing was performed for within-run (using three donors), within-device (20-day tests using QC Levels 1 and 2), between lots (using three different cartridge lots), and between-operator (using five operators). Also, multi-site reproducibility tests were performed across three sites, three lots, three operators and three analyzers. Results included Coefficient of Variance (CV) values for all precision tests. The CV of the test results for the Functional Fibrinogen Assay was less than 10% for the MA and FLEV parameters.

CFF Linearity

The Functional Fibrinogen assay was tested for linearity according to CLSI EP6-A. The MA parameter was found to be linear with respect to concentrations of Fibrinogen, as shown below.

Packaging Symbols

Symbol	Description
	Manufacturer (address for)
	Authorized representative in the European community (address for)
	CE mark
	Catalogue number
	Batch code
	Control
	In vitro diagnostic device
	Consult instructions for use
	Do not reuse
	Temperature limitation for storage
	(Applies to USA only) USA Federal Law restricts the sale, distribution or use of this device to, by or on the order of a physician.
	Contains sufficient for <n> tests
	Use by <date>

Symbol	Description
	Do not use if package is damaged
	Type of venous blood collection tube to be used for the test. Color-coded top conforms to ISO 6710:1995

Legal Notice

Licensed for sale or use only in in-patient hospitals and any laboratories that are primary providers of such hospitals.

Bibliography

1. CLSI EP5-A2 Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition
2. CLSI EP6-A Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline
3. CLSI EP7-A2 Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition
4. CLSI C28-A3c Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition

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