# CoaguChek<sup>®</sup> XS

## User's Manual Manuel d'utilisation

On the packaging and on the type plate of the meter you may encounter the following symbols and abbreviations shown here with their meaning:

Use by

Batch code/ Lot number

In vitro diagnostic medical device

This product fulfills the requirements of the European Directive 98/79/EC on in vitro diagnostic medical devices.

Catalogue number

Consult instructions for use

Caution, consult accompanying documents. Refer to safety-related notes in the manual accompanying this instrument.

Manufacturer

Temperature limitation (Store at)

The system fulfills the Canadian and U.S. safety requirements (UL LISTED, in accordance with UL 61010A-1:02 and CAN/CSA-C22.2 No. 61010-1-04)

## The CoaguChek XS System

The CoaguChek XS System is used for quantitative monitoring of coagulation values (common terms: thromboplastin time, PT, Quick value) with CoaguChek XS PT test strips.

Self-monitoring and self-adjustment of oral coagulation therapy using coagulation values determined with the aid of the CoaguChek XS System may only be undertaken after consultation with your doctor and comprehensive instruction by a qualified healthcare professional. Your pharmacist / specialist supplier can provide you the training and/or put you in touch with associations or institutions offering training in coagulation self-monitoring. Please read the special notes on INR self-monitoring in the package insert accompanying the test strips.

### WARNING

---

The CoaguChek XS Softclix lancing device is intended only for patient self-monitoring by a single person. **It must not be used to collect blood in a multi-patient setting as it does not incorporate any features to guard against cross-infection.**

---

This manual contains all the information you need to operate and care for the CoaguChek XS System. Please **read this entire manual carefully** before you use the meter.

**Healthcare professionals:** In addition, observe carefully the precautions and procedures applying to professional use provided in the relevant sections of this manual.

## Revision History

<b>Manual version</b>	<b>Revision date</b>	<b>Changes</b>
1.0	2009-09	New document
2.0	2010-03	Cleaning/disinfection recommendations and presentation safety messages revised. Screens updated.
3.0	2011-04	Change in lancing device. Memory: 300 test results with date and time. Screens updated.
4.0	2015-05	Manual update: updated cleaning/disinfection section, minor revisions

<b>The CoaguChek XS System</b>	3
<b>Introduction</b>	9
The CoaguChek XS System.....	9
Test principle.....	10
Contents of the pack.....	10
Operating conditions.....	11
Quality control.....	12
<b>Important Safety Instructions and Additional Information</b>	13
<b>The CoaguChek XS Meter</b>	14
The CoaguChek XS meter.....	15
Batteries.....	16
<b>Starting Up</b>	17
Inserting the batteries.....	18
<b>Code Chip</b>	21
Inserting the code chip.....	22
<b>Meter Setup</b>	24
Settings summary.....	25
Setting up the meter (Set Mode).....	26
Setting the date format.....	29
Setting the date.....	30
Setting the time format.....	32
Setting the time.....	33
Selecting the unit.....	34
Setting the beep tone.....	35
Setting the therapeutic range (INR).....	36

<b>Testing a Capillary Blood Sample</b>	41
Important notes.....	42
Getting a good capillary blood sample .....	44
Preparing to test.....	45
Performing a test.....	46
Notes concerning the results display .....	57
<b>Testing by Healthcare Professionals</b>	61
Getting a good capillary blood sample in a healthcare setting	62
Disposal by healthcare professionals .....	62
<b>Memory</b>	63
Viewing test results.....	63
Erasing memory .....	66
<b>Cleaning and Disinfecting the Meter</b>	69
What is the difference between cleaning and disinfecting? ...	70
When should the meter be cleaned and disinfected? .....	70
Recommended cleaning/disinfection agents.....	71
Cleaning/disinfecting the meter housing .....	72
Cleaning/disinfecting the test strip guide .....	73
<b>Error Messages</b>	77
Error messages overview .....	78
Error messages after turning the meter on .....	80
Error messages when preparing to test .....	82
Error messages after the code number was confirmed .....	86
Error messages during or after blood application .....	88
Error message during infrared downloads .....	96
<b>Display Information and Symbols</b>	98

<b>Further Information</b>	103
Ordering .....	103
Product limitations .....	103
<b>Product Specifications</b>	104
Operating conditions and technical data .....	104
Sample material .....	105
Storage and transport conditions .....	105
Disposal of the CoaguChek XS meter .....	106
Disposal of used batteries .....	106
Information service .....	107
Repairs .....	107
<b>Guarantee</b>	108
<b>Alphabetical Index</b>	109

This page intentionally left blank.

# Introduction

## The CoaguChek XS System

The CoaguChek XS System (CoaguChek XS meter and CoaguChek XS PT test strips) quantitatively determines prothrombin time (PT/Quick value/ INR) using capillary blood from a fingertip or untreated venous whole blood. The CoaguChek XS System makes coagulation testing easy. The CoaguChek XS meter guides you through the test step by step using the symbols in the display. The code chip accompanying the test strips contains lot-specific information about these test strips, including calibration data that allow the correct results to be calculated, and the expiry date for the test strips. You only need to insert the code chip, turn the meter on, insert the test strip, and apply a blood sample. When the CoaguChek XS meter is ready for a measurement, the result will be displayed about 1 minute after application of the sample. After the measurement, the meter automatically stores the result in memory.

If you have any questions about the CoaguChek XS meter, please contact your local customer support and service centre. You will find contact details on page 107.

**Please note:** Before you use the meter for the first time (i.e. after you have first inserted the batteries), you must set the date and time correctly to allow you to carry out measurements properly. Each time you replace the batteries you need to check (and, if necessary adjust) the date and time.

### **Test principle**

The CoaguChek XS PT test strip contains a lyophilized reagent (reagent in dried form). The reactive components of this reagent consist of thromboplastin and a peptide substrate. When a sample is applied, thromboplastin activates coagulation, which leads to the formation of thrombin. At the same time the meter starts to measure the time. The enzyme thrombin cleaves the peptide substrate, generating an electrochemical signal. Depending on the time elapsed when it first appears, this signal is then converted by means of an algorithm into customary coagulation units (INR, %Quick, seconds) and the result is displayed.

### **Contents of the pack**

- CoaguChek XS meter
- 4 alkali-manganese batteries, 1.5 V, type AAA (LR03)
- CoaguChek XS Softclix lancing device with instructions for use
- Lancets for the CoaguChek XS Softclix lancing device
- System carry-case
- User's Manual

## Operating conditions

To ensure that your CoaguChek XS System functions properly, please observe the following guidelines:

- Only use the meter at a room temperature between 15 °C and 32 °C.
- Only use the meter at a relative humidity between 10% and 85%.
- When testing, place the meter on a level, vibration-free surface or hold it so it is roughly horizontal.
- If the meter is to remain unused for a longer period of time, keep it in the carry-case supplied.
- Use the meter only up to a maximum altitude of 4300 metres.

---

### **Electromagnetic interference**

Strong electromagnetic fields may interfere with the proper operation of the meter. Do not use the meter near strong electromagnetic fields.

---

### **Quality control**

The CoaguChek XS System has a number of inbuilt quality-control functions such as:

- A check of the electronic components and functions every time the meter is turned on.
- A check of the test strip temperature while a test is in progress.
- A check of the expiry date and lot information on the test strip.
- A quality control function is incorporated into the test strip. Quality control and system checks using control test solutions that you may be familiar with from other systems are no longer required.

## Important Safety Instructions and Additional Information

This section explains how safety-related messages and information related to the proper handling of the system are presented in the CoaguChek XS User's Manual. Read these passages carefully.

---

The safety alert symbol by itself (without a signal word) is used to promote awareness to hazards which are generic or to direct the reader to related safety information

---

**WARNING**

---

Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

---

**CAUTION**

---

Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

---

**NOTICE**

---

Indicates a hazardous situation which, if not avoided, may result in damage to the system.

---

- Important information that is not safety relevant is presented against a colored background (without a symbol). Here you will find additional information on correct use of the meter or useful tips.
- Instructions that require you to actively do something also appear on a blue background.

## The CoaguChek XS Meter

## The CoaguChek XS meter

- A Display**  
Shows results, information, symbols and results recalled from memory.
- B M (memory) button**  
Press this button to recall results from memory and to change the meter settings. You also press this button to confirm the code number displayed before each test.
- C On/Off button**  
Press this button to turn the meter on or off.
- D Test strip guide cover**  
Remove this cover to clean the test strip guide.
- E Test strip guide**  
Insert the test strip here.
- F Battery compartment cover**  
Covers the battery compartment (4 1.5 V alkali-manganese batteries type AAA (LR03)).
- G Code chip slot**  
Insert the code chip here.
- H Set button**  
The Set button is located on the left side of the meter. Press this button to enter or change the meter settings.
- I Infrared window**  
You can transfer data stored in memory via this infrared interface.

### Batteries

To save power, the CoaguChek XS meter automatically turns itself off after 3 minutes unless a button has been pressed or a new test strip has been inserted. When the meter turns itself off, all results obtained up to that point remain in memory. The meter briefly indicates the battery power level when it is turned on. The battery symbol is divided into four segments which correspond to the battery power level.

When replacing the batteries you must insert the new batteries within one minute of removing the old ones, to keep the date and time settings. If you take longer than this, you must re-enter the date and time. Use only alkali-manganese batteries type AAA (LR03).

The CoaguChek XS meter is shipped with high quality batteries. Since quality varies among different battery brands, Roche recommends replacing your used batteries with new batteries of a similar high quality and not mixing different brands. Do not mix new and already used batteries.

Please be aware that battery life may be affected by many factors, such as battery quality, operating conditions (e.g. ambient temperature), frequency of use, and test duration. Therefore replace used batteries after one year at the latest.

Results are retained in memory together with the date and time even when no batteries are inserted. All other settings are retained, too.

Think of the environment. Dispose of used batteries sensibly.

**WARNING**

---

Do not throw batteries onto a fire. They may explode!

---

## Starting Up

Before using the meter for the first time, carry out the following steps:


- 1 Insert the batteries
- 2 Insert the code chip (you can also carry out this step immediately before testing)
- 3 Set the current date and time
- 4 Select the coagulation unit in which the results are to be displayed

### Inserting the batteries

- 1** With the meter turned off, turn it over.
- 2** Gently press the tab on the battery compartment cover towards the centre of the meter and lift off the cover.
- 3** Insert the four batteries in the battery compartment as indicated. Pay attention to the positions of “+” (top of battery) and “-” (flat end). Once you have inserted the batteries, the meter turns itself on after about 5 seconds.

Only use alkaline manganese batteries (1.5 V, AAA, LR03). Think of the environment. Dispose of used batteries sensibly.



- 4 Close the battery compartment.
- 5 Check that all of the display segments and symbols are displayed properly. If your display is not functioning properly, you must not perform further tests as results may be misread if a segment is missing. Contact your local customer support and service centre in this case.

### Tip

The meter performs an automatic display check when it is turned on. If the display check is too quick, you can put it “on hold” by pressing and holding down the On/Off button (ⓘ) the next time you turn the meter on.

For as long as you keep the button pressed, the display remains frozen.

- 6** On the next screen that appears check that the date and time are correctly set.

If these have not yet been set or if they have been lost (because the batteries were removed from the meter for more than one minute), your CoaguChek XS meter now automatically enters Set Mode. After you have set the date and time, the meter enters Test Mode. If you wish to enter any further settings, turn to the chapter *Setting up the meter (Set Mode)* starting on page 26.

- 7** Turn the meter off again when the settings are correctly entered.

## Code Chip

The code chip provides the meter with important information which it needs to perform the coagulation test. The chip contains information about the test method, the lot number and the expiry date. The meter is ready to use once the code chip has been inserted.

- Each time you test, check that the correct code chip is inserted in the meter.
- When you start each new pack of test strips you must remember to remove the old code chip from the meter and replace it with the new code chip supplied with the test strips.
- Each code chip belongs to a particular lot of test strips. Only remove the code chip when you are testing with test strips taken from a new pack (with a new code chip).
- Protect the code chip from moisture and equipment which produces magnetic fields.

### Inserting the code chip

- 1 Remove the old code chip, if one is inserted in the meter.

Discard the old code chip with your household waste.

- 2 Always make sure that the number on the code chip matches the number on the label of the test strip container.

**CAUTION**

---

Using the wrong code chip can produce incorrect results.

---

- 3 Slide the new code chip into the slot on the side of the meter as shown until you feel it snap into place.

Every time you insert a test strip in the meter, the display shows the number of the code chip that is presently inserted. Always compare the code number you see on the display with the number that is printed on the test strip container. If the two numbers are identical, confirm by pressing the **M** button (see page 49).

If the code numbers are not identical, turn the meter off and insert the correct code chip. Discard the old code chip to prevent mix-ups.

If the code chip is missing or incorrectly inserted, **Error** and **Code** appear in the display (please refer also to the chapter *Error Messages* starting on page 77).

## Meter Setup

**Please note:** If you have **not** set the date (after turning on for the first time or because the batteries were removed from the meter for more than one minute), you cannot perform a test. In that case turning on the meter takes you immediately to Set Mode, where you must set the date and time. Once you have set date and time the meter automatically enters Test Mode.


**Settings summary**

Setting	Options	Presetting *
Date format	Day-Month-Year (31-12-00) Month-Day-Year (12-31-00) Year-Month-Day (00-12-31)	31-12-00
Date		31-12-00
Time format	24-hour format (24h) 12-hour format (12h) with a.m./p.m.	24h
Time		12:00
Unit	%Q Sec INR	INR
Beep tone	On Off	On
Therapeutic range (only INR)	On Off	Off
Set range [INR]	Lower range (1.5 – 3.5) Upper range (2.5 – 4.5)	1.5 INR 2.5 INR

\* “Presetting” means the setting that is in the meter when it leaves the factory.

## Setting up the meter (Set Mode)


All settings are entered with the **Set** button and the **M** button.

- 1** If the meter did not automatically enter Set Mode when you turned it on, press the **Set** button (on the left side of the meter) to put it in Set Mode.

You can also press the **Set** button to change to Set Mode when the meter is already turned on.

- 2** If the currently displayed setting is correct (e.g. the date is correct and you only wish to change the time), you can press the **Set** button to move on to the next setting, **or**

- 3** Press the **M** button to change the setting that is currently flashing. You may press the **M** button repeatedly (or press it and hold it down) until you have reached the correct setting. Settings where there are only two options (time format, beep tone) can be toggled on and off with the **M** button.
- 4** Press the **Set** button again to confirm the current setting (which is saved to memory) and move on to the next one.

- 5** You can only move forward from one setting to the next, not backward. Corrections can only be made by repeating the settings. You can exit the settings procedure at any time by pressing the On/Off button . All settings you have already entered are retained.

## Setting the date format


First, select the date format (**the entire date display flashes**).

The display options are as follows:

- **31-12-00 (= presetting) Day-Month-Year**
- 12-31-00 Month-Day-Year
- 00-12-31 Year-Month-Day

When elements in the display are surrounded by a halo, it indicates that they are flashing in the meter's display.



The **Set** button is located on the left side of the meter. Pressing the **Set** button takes you to Set Mode. This is indicated by the word "**set**" in the display.


- 1** Press the **M** button to select the date format. 31-12-00 flashes in the display. You can now select the date format.
- 2** Press the **Set** button to save this setting and continue with the date setting. The display automatically moves on to the date setting.

## Setting the date

The next three settings involve entering first the **year**, then the **month**, and finally the **day**.


- 3 Press the **M** button to change the number that is currently flashing and set the current year.
- 4 Press the **Set** button to save the year. The display then automatically moves on to the month setting.
- 5 The preset month flashes. Press the **M** button until the correct month is displayed.
- 6 Press the **Set** button to save the setting. The display then automatically moves on to the day setting.

- 7** The preset day flashes. Press the **M** button until the correct day is displayed.
- 8** Press the **Set** button to save the setting. The display then automatically moves on to the time format setting.

## Setting the time format

Now select the time format. You can choose between the 24-hour time format (presetting) and the 12-hour time format with “a.m.” or “p.m.”.



- 9** Press the **M** button to toggle between 24-hour and 12-hour format.
- 10** Press the **Set** button to save the correct setting and continue with the time setting. The display automatically moves on to the time setting.

## Setting the time

First you set the hour, then the minutes.

- 11 Press the **M** button to change the number that is currently flashing. Press the **Set** button to confirm. You can now set the minutes (again using the **M** button).

If you have selected the **12-hour** time format, “**a.m.**” changes to “**p.m.**” and vice versa when you reach “12:xx”.

- 12 Press the **Set** button to save this setting and continue with selection of the unit for reporting the result.

## Selecting the unit

Choose the coagulation unit you wish the meter to use to display the result. The options are:


- %Q
- Sec
- **INR (presetting)**

### WARNING

Ask your doctor if you are not sure which unit to select.

The current coagulation unit flashes.

**13** Press the **M** button to cycle through the three options.

**14** Press the **Set** button to save the correct unit. The display then automatically moves on to the beep tone setting.


If you have selected “INR” as your coagulation unit, you have the option of entering your therapeutic target range. If you enter this, you will be able to see after every test whether the result is within your therapeutic target range (see page 36).

## Setting the beep tone

After you have set the unit, you can choose whether you want the beep tone “On” or “OFF”. If you have selected “On,” the meter will beep in the following situations:

- when it detects a test strip,
- when pre-heating of the test strip is complete and you need to apply a sample,
- when it detects a sample,
- when the result is displayed and
- if an error occurs (three short beeps).

We recommend you leave the beep tone turned on.


- 15** Press the **M** button to toggle between “OFF” and “**On**” (“On” is the presetting).
- 16** Press the **Set** button to save the chosen setting. The display automatically moves on to the next setting option.

## Setting the therapeutic range (INR)

What happens next depends on the coagulation unit you have chosen for your results. If you selected “%Q” or “Sec”, the settings are now complete, and “End” appears in the display (see page 38).



If you selected “INR” as the unit, you can now choose whether you want a further element to be displayed if the result is outside (above or below) your therapeutic target range. If you switch on this display option, an up or down arrow is displayed whenever the result is above or below the target range you have set. This highlights the fact that the result is outside your target range.

**17** Press the **M** button to toggle between “**OFF**” and “On”.

**18** Press the **Set** button to save the chosen setting and continue with the meter setup.

If you have activated this display option (“On”), you are now able to enter the limits of your target range. Consult your doctor about what target range is appropriate for you.

**19** Press the **M** button to set the **lower** limit (symbolized by the down arrow) within the range 1.5 – 3.5 INR.


**20** Press the **Set** button to save the chosen lower limit and continue with the upper limit setting.

**21** Now press the **M** button to set the **upper** limit (symbolized by the up arrow) within the range 2.5 – 4.5 INR. The selectable range will be at least 0.1 INR above the lower limit you have selected.


**22** Press the **Set** button to save the chosen upper limit and close the meter setup.

## Meter Setup

The meter setup is now complete.

This display screen appears automatically after you have completed the last setting and remains for a few seconds to indicate that the setup procedure is now complete.

After that the meter automatically enters Test Mode.

If you do not wish to perform a test at this point, turn the meter off.

You can exit the meter setup procedure at any time (e.g. if you have made an incorrect entry) by pressing the On/Off button . Then enter the meter setup again by turning on the meter with the **Set** button, and press the **Set** button repeatedly until you reach the setting you wish to change.

## Meter Setup

This page intentionally left blank.

## Testing a Capillary Blood Sample

### What you need

- CoaguChek XS meter
- The code chip already inserted in the meter (each test strip container has its own code chip)
- Test strips belonging to the code chip mentioned above
- Lancing device (e.g. CoaguChek XS Softclix)
- Lancet (e.g. CoaguChek XS Softclix Lancet)
- Lint-free cloth or paper tissue

### **Important notes**

#### **Always ...**

- ... close the container immediately after having removed a test strip.
- ... operate the meter at an ambient temperature between 15 °C and 32 °C.
- ... place the meter on a level, stable surface (table) or hold it so it is roughly horizontal.
- ... follow the information on correct handling of test strips in the package insert.
- ... keep the test strip guide and housing clean. Refer to the chapter *Cleaning and Disinfecting the Meter*.

**Never ...**

- ... store the meter at extreme temperatures.
- ... store the meter in damp or humid conditions without protection.
- ... remove or insert the code chip while the meter is performing a test.
- ... use the code chip from a pack other than the one in use.
- ... touch or remove the test strip during a test.
- ... wait more than 15 seconds after lancing the fingertip before applying the blood.
- ... add more blood after the test has begun.
- ... perform a test with a drop of blood from a previous puncture.

---

Failure to comply with the above may lead to inaccurate results. An incorrect result may lead to an error in diagnosis, therefore posing danger to the patient.

---

## Getting a good capillary blood sample

To get a suitable drop of blood:

- Wash your hands in warm water and dry thoroughly.

### WARNING

---

Residues of water on the skin can dilute the drop of blood and so produce false results.

---

- Let your hand hang at your side before lancing your finger.
- Immediately after lancing, massage gently along the side of your finger to obtain a sufficiently large blood drop without pressing or squeezing too hard.

### WARNING

---

The CoaguChek XS Softclix lancing device is intended only for patient self-monitoring by a single person. **It must not be used to collect blood in a multi-patient setting as it does not incorporate any features to guard against cross-infection.**

---

## Preparing to test

- 1 Have the test strip container to hand.
- 2 Make sure that the code chip belonging to these test strips is inserted in the meter.
- 3 Prepare the lancing device by inserting a fresh lancet.

Do not lance your finger until asked to do so further on in this description.

## Performing a test

- 1 Wash your hands with soap and warm water. Dry them thoroughly.



### WARNING

Residues of water on the skin can dilute the drop of blood and so produce false results.

- 2 Place the meter on a level, vibration-free surface or **hold it in your hand so it is roughly horizontal**. Turn the meter on by pressing the On/Off button (M). Alternatively, you can insert a test strip to turn it on.
- 3 Check that all of the display symbols are displayed properly as shown on the picture. Results may be misread if a segment is missing.


- 4 Check the battery level. If there are no bars left in the battery symbol, you cannot perform any more tests.
- 5 Check that the date and time are correct. Correct any wrong entries as described from page 26 onwards.
- 6 The flashing test strip symbol prompts you to insert a test strip. Remove a test strip from its container. **Immediately after removing a test strip, close the container again with the stopper.**

**NOTICE**

Exposure to external influences (e.g. humidity) may deteriorate the test strips and may lead to error messages.


## Testing a Capillary Blood Sample

- 7** Hold the test strip so the lettering, the arrows and the blood drop symbols are facing upward.
- 8** Slide the test strip into the test strip guide in the direction indicated by the arrows.

Slide the test strip in as far as it will go. A beep tone indicates that the meter has detected the test strip (provided the beep tone is turned on in the settings).

The code number of the code chip inserted in the meter flashes in the display. Make sure that this number is identical with the code number printed on the test strip container.

**9** If the two numbers are identical, confirm by pressing the **M** button.

The code number stops flashing.


If the code numbers are not identical, remove the wrong code chip and insert the code chip that was supplied with the test strips you are using.


**CAUTION**

Using the wrong code chip can produce inaccurate results.

## Testing a Capillary Blood Sample

The hourglass symbol shows that the test strip is warming up. When the warming-up process is complete, a further beep tone (provided the beep tone is turned on) indicates that you can now apply blood.

The blood drop symbol and the application area flash to indicate that the meter is ready to perform the test and is waiting for blood to be applied.

At the same time a 180-second countdown begins. You must apply the drop of blood to the test strip within this time, otherwise you will receive an error message (**error**, and **"000"** where the result would normally be displayed).

**10** Now lance the side of a fingertip with the lancing device.

We recommend obtaining the capillary blood from the side of the fingertip as this causes least pain.

Massage the lanced finger until a drop of blood is formed.

  
**CAUTION**

---

**Do not press or squeeze the finger.**

---

  
**WARNING**


**Apply the first drop of blood from the finger.**

---

The CoaguChek XS Softclix lancing device is intended only for patient self-monitoring by a single person. **It must not be used to collect blood in a multi-patient setting as it does not incorporate any features to guard against cross-infection.**

---

## Testing a Capillary Blood Sample

- 11** Apply the blood directly from the finger to the semicircular, transparent sample application area of the test strip.

Alternatively, you can touch the blood drop against the side of the sample application area (instead of applying it to the centre from above). The test strip draws up the blood by capillary action.

During this process you must hold the blood drop to the test strip until the flashing blood drop symbol has disappeared and the meter beeps (provided the beep tone is turned on).


---

Apply the blood drop to the test strip **within 15 seconds** of lancing the fingertip. Applying blood after this period of time would falsify the result, as the coagulation process would already have begun.


---

## Testing a Capillary Blood Sample


You hear a beep tone when you have applied enough blood (provided the beep tone is turned on). The blood drop symbol disappears and the test starts.

Do not add more blood. Do not touch the test strip until the result is displayed.


The meter now performs an automatic quality control check on the test strip. "QC" appears in the display.

Following a successful outcome of the quality control check, a tick appears after "QC".

It is at this point that the coagulation measurement begins. The hourglass flashes until the meter has the result.

The result is displayed in the unit you chose when setting up the meter. It is automatically saved to memory. You can change the unit if you wish by pressing the **Set** button while the result is displayed (see page 58). The **Set** button is located on the left side of the meter.

Enter the test result in your patient logbook. When interpreting your result, please also read the test strip package insert carefully.

- 12** Remove the test strip from the test strip guide.
- 13** Turn the meter off.
- 14** Discard the used test strip with your regular household waste.

### **WARNING**



---

Dispose of used lancets carefully (e.g. use a sturdy sharps container with lid) so that the needles cannot cause injury to yourself or to others.

---


- 15** Clean the meter if this becomes necessary (see page 69).

**Notes concerning the results display**

If your results are displayed in INR and you have activated the relevant option, values that are outside a defined target range are flagged with an up arrow (too high) or a down arrow (too low).

## Testing a Capillary Blood Sample

A “c” is displayed along with the result.


This may occur if the haematocrit value is very low or due to erroneous blood collection (e.g. wet hands). Repeat the measurement and make sure that your hands are dry. If the message persists, please contact your healthcare professional for a haematocrit check.

You can temporarily change the unit in which the result is presented by pressing the **Set** button while the result is displayed.

Each press of the button takes you to the next available unit (%Q - Sec - INR - ...).

When you next perform a test, the result is again displayed in the standard unit you selected during setup.

The valid measuring ranges for CoaguChek XS PT test strips are:

- INR: 0.8–8.0
- %Q: 120–5
- Sec: 9.6–96

Results that are outside the measuring range are indicated by the symbols > (greater than) or < (smaller than) and flashing. In this case, refer to the section *Test limitations and known interferences* in the package insert supplied with the test strips, and repeat the test.

## Testing a Capillary Blood Sample

This page intentionally left blank.

## Testing by Healthcare Professionals

### Protection against infection

There is a potential risk of infection. Healthcare professionals using the CoaguChek XS system must be aware that any object coming into contact with human blood is a potential source of infection. Healthcare professionals must also be aware that any cross-contamination is a potential source of infection for patients. (See: Clinical and Laboratory Standards Institute: Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline – Third Edition; CLSI document M29-A3, 2005).

- Use gloves.
- Use a separate lancet / separate lancing device for each patient.
- Dispose of used lancets in a sturdy sharps container with lid.
- Dispose of used test strips according to your institution's infection control policy.
- Follow all health and safety regulations in force locally.
- Use only a lancing device that is approved for use by healthcare professionals (for example, the Accu-Chek Safe-T-Pro Plus lancing device from Roche shown below). Follow the manufacturer's instructions for use in a multi-patient setting.

## Getting a good capillary blood sample in a healthcare setting

---

### CAUTION

When washing or disinfecting the patient's finger, allow it to dry thoroughly. Residues of water or disinfectant on the skin can dilute the drop of blood and so produce false results.

---

- Select and lance a puncture site along the side of a fingertip.
- Use gentle pressure, massage the finger towards the fingertip to encourage a drop of blood to form (at least 8  $\mu$ L).
- Follow the procedure for testing as outlined from page 45.
- Use only a lancing device that is approved for use by healthcare professionals (for example, the Accu-Chek Safe-T-Pro Plus lancing device from Roche).

## Disposal by healthcare professionals

Dispose of used items in line with the disposal policy of your hospital, institute or medical practice.


## Memory

The CoaguChek XS meter has a 300-value memory for saving results together with time and date.

### Viewing test results


- 1 Place the meter on a level, vibration-free surface or hold it in your hand so it is roughly horizontal.
- 2 Press the **M** button to turn the meter on in Memory Mode.

The meter runs the usual display check and displays the most recent result. The display also shows **mem**, signifying that this is a result from memory and not one that has just been obtained by testing.


The date and time of testing are also shown.

- 3 Press the **M** button again to move one result back in time.


When you have reached the oldest result, the meter displays a blank entry.

When all memory locations are full, the oldest result is erased to make way for the most recent one.

If no result has yet been saved, the display looks similar, but has a "0" in the top right corner.

If you press the **Set** button while viewing a result in memory, you can temporarily change the unit in which the result is displayed (in the same way as directly following a test). The **Set** button is located on the left side of the meter. Each press of the button takes you to the next available unit (see page 58).

Results can be downloaded from memory via an inbuilt infrared interface.

**Please note:** When downloading results via the infrared interface, avoid strong stray light, as this can cause the download to break off (see page 96).


For more information on downloading results to a suitable system (such as a PC), please call your local customer support and service centre (see page 107).

### Erasing memory

You can erase the entire memory from your CoaguChek XS meter. You cannot, however, erase individual results.

- 1** Place the meter (which must be turned off) on a level, vibration-free surface or hold it so it is roughly horizontal.
- 2** To turn the meter on,
  - first press **and hold down** the **M** button,
  - then press the **On/Off** button as well, and
  - hold down both buttons for at least 5 seconds.


The meter runs the usual display check and then displays **mem** (flashing) and **clr**. The number of results in memory is shown in the top right corner of the display.


- 3 Press the **M** button to confirm that you wish to erase the entire memory.

To exit Erase Mode without erasing the results, press the **On/Off** button instead.

## Memory

The hourglass flashes while the results are being erased.

The counter is set to "0" and three bars are displayed, indicating that there are no results saved in memory (the entire memory has been erased).

## Cleaning and Disinfecting the Meter

---

**NOTICE** To avoid malfunction of the meter, follow the procedures below to clean and disinfect the meter.

---

- Make sure the meter is turned off.
- Do not use sprays of any sort.
- Ensure that the cloth or cotton swab/bud is only damp, not wet, to protect against moisture entering the meter.

**Useful tip for self-testers:** To prevent soiling, apply blood via side-dosing directly from fingertip.

### **What is the difference between cleaning and disinfecting?**

Cleaning is the removal of dirt from the meter.

Disinfecting is the removal of most, but not all, disease-causing and other types of microorganisms (bloodborne pathogens) from the meter.

### **When should the meter be cleaned and disinfected?**

- Clean and disinfect the meter to remove visible dirt or other material on a regular basis.
- Clean and disinfect the meter always when there is blood on it.
- Clean and disinfect the meter before allowing anyone else to handle the meter, for instance, if you have someone assisting you.
- If the meter is used in a professional environment, the exterior of the meter and the test strip guide should be cleaned and disinfected before being used on another patient.

## Recommended cleaning/disinfection agents

---

### Healthcare professionals:

- Observe the disinfection guidelines of your institution.
  - Use gloves.
- 

Gently wipe the meter's surface with one of these cleaning agents:

- A soft cloth slightly dampened (not wet) with a small amount of liquid soap diluted in water
- Rubbing alcohol (70% ethanol or isopropyl alcohol)
- Alcohol-based disinfectant (a mixture of 1-propanol (400 mg/g), 2-propanol (200 mg/g) and glutaraldehyde (1.0 mg/g))
- 10% sodium hypochlorite solution (1 part bleach to 9 parts de-ionized water, made fresh every 24 hours)
- Disposable wipes (containing quaternary ammonium compounds up to 0.5 % (single compound or mixture) in isopropyl alcohol (isopropanol) up to 55 %)

Ask your healthcare professional or pharmacist to recommend commercially available wipes containing the acceptable cleaning agents.

### Cleaning/disinfecting the meter housing

Use the solutions recommended on page 71 for cleaning/disinfecting the meter exterior. Apply the solutions for a contact time of > 1 minute (refer to the corresponding product labeling).

- Ensure that the blue test strip guide cover remains tightly closed while cleaning the meter housing.
- Do not let liquid accumulate near any opening.

#### NOTICE

---

Ensure that **no** liquid enters the meter when you clean the exterior. If moisture enters the meter, it may cause malfunction of the instrument.

---

- 1 With the meter turned off, wipe the meter's exterior clean.
- 2 With a fresh dry cloth or lint-free tissue, wipe away residual moisture and fluids after cleaning the housing.
- 3 Allow wiped areas to **dry for at least 10 minutes** before performing a test.

### **Cleaning/disinfecting the test strip guide**

Use the solutions recommended on page 71 for cleaning/disinfecting the test strip guide. Apply the solutions for a contact time of > 1 minute (refer to the corresponding product labeling) using lint-free cotton swabs/buds.

- 1** Remove the test strip guide cover to clean it. (Use your thumbnail to open the cover of the test strip guide by pressing its front edge upward.) Move the cover safely away from the meter.
- 2** Rinse the cover with warm water or wipe it clean using the solutions recommended above. Let the test strip guide cover dry for at least 10 minutes before re-attaching it.

## Cleaning and Disinfecting the Meter

- 3** Hold the meter upright with the test strip guide facing down.
  - Clean the easily accessible white areas with a moistened cotton swab/bud.
  - Ensure the swab/bud is only damp, not wet. Wipe away residual moisture and fluids.

**NOTICE**

---

**Damage to the instrument**

- Ensure that **no** liquid enters the meter. If moisture enters the meter, it may cause malfunction of the instrument.
  - Do not insert any objects in the test strip guide. Doing so may damage the electrical contacts behind the test strip guide.
- 

On completion of cleaning:

- 4 With the cover off, let the test strip guide dry for **at least 10 minutes**.
- 5 After this time, re-attach the test strip guide cover to the housing. Make sure that the cover is properly closed. You will hear it snap into place.

## Cleaning and Disinfecting the Meter


This page intentionally left blank.

## Error Messages

Depending on the circumstances an error message may appear on the display of your meter. Your first step on receiving an error message should be to try the solution suggested for that error. If the problem persists, please contact your local customer support and service centre (see page 107).

A complete list of display messages and symbols is provided in this User's Manual starting on page 98.



## Error messages overview

Error message (always prefaced by <b>error</b> )		Brief description	Refer to page
Test strip symbol flashing		Test strip in meter or test strip unusable or the test strip is not a CoaguChek XS PT test strip	80 84
Temperature symbol flashing		Meter too warm or too cold	80
Battery symbol flashing		Battery level too low	82
Cover symbol flashing		Test strip guide cover open	82
IR symbol flashing		Infrared download error	96
QC symbol flashing	<b>QC</b>	Quality control: test strip unusable	90
Code symbol flashing	<b>code</b>	Code chip error	84

Error message (always prefaced by <b>error</b> )	Brief description	Refer to page
Error number: <i>000</i>	Allowed time for blood application exceeded	88
Error number: <i>3</i>	Beyond expiry date	86
Error number: <i>4</i>	Test strip unusable or internal error	86
Error number: <i>5</i>	Blood application error	88
Error number: <i>6</i>	Measurement error	90
Error number: <i>7</i>	Measurement error	92
Error number: <i>8</i>	Internal error	94
Error number: <i>9</i>	Internal error	94

Error numbers 1 and 2 are reserved for internal purposes by Roche.

## Error messages after turning the meter on

Display	Description of error
 <p>The display shows the time 11:23, the date 15-05-15, and the word 'error'. Below the display is a battery icon with a full charge indicator.</p>	<p><b>Error: test strip</b> A test strip was already inserted when the meter was turned on with the <b>On/Off</b> button.</p>
 <p>The display shows the time 11:23, the date 15-05-15, and the word 'error'. Below the display is a thermometer icon.</p>	<p><b>Error: ambient temperature</b> The meter is too cold or too warm to measure correctly.</p>


**Solution**

Remove the test strip. The error message disappears, the meter is ready.

Turn the meter off and allow it to stand for about 30 minutes at room temperature (+15 °C to +32 °C).

**Display**


**Description of error**

**Error: battery**

Battery level too low

**Error messages when preparing to test**

**Error: cover**

The test strip guide cover is not properly closed.


**Solution**

Insert new batteries as described on page 18.

Close the test strip guide cover. The error message disappears, the meter is ready.


**Display**

**Description of error**

**Error: test strip**

The test strip is unusable or the test strip is not a CoaguChek XS PT test strip.

**Error: code chip**



The code chip is missing, not properly inserted in the meter, or damaged.

**Solution**

Remove the test strip and repeat the test with a new or correct test strip.

Insert the code chip, or remove it and re-insert it.

## Error messages after the code number was confirmed

Display	Description of error
 A digital display showing the number 3. The number is formed by seven segments. Above the number, the word "error" is displayed. The display is surrounded by a square frame with tick marks at the corners and midpoints of the sides.	<p><b>Error: test strip</b> The test strip is beyond its expiry date.</p>
 A digital display showing the number 4. The number is formed by seven segments. Above the number, the word "error" is displayed. The display is surrounded by a square frame with tick marks at the corners and midpoints of the sides.	<p><b>Error: test strip or meter</b> The test strip is unusable, or the meter is defective.</p>

### Solution

First check whether the date setting is correct in the meter. If it is not, set the correct date.


If the date is correct: turn the meter off, remove the code chip and test strip. Use a test strip from a new lot of test strips and insert the code chip that came with the new lot.

Turn the meter off, remove the test strip and re-insert it. If the error message appears again, discard the unusable test strip and use a new one.

If you still get the error message after you have inserted a new test strip, turn the meter off again. Wait two minutes before turning the meter on. If the error message still appears, the meter might be defective. Contact your local Roche Diagnostics customer support and service centre in this case (see page 107).

**Display**


**Description of error**

**Error: time exceeded**

The maximum allowed time of 180 seconds for applying blood has been exceeded.

**Error messages during or after blood application**

**Error: blood application**



Error applying blood to the test strip.

**Solution**

Turn the meter off, remove the test strip and start a new test using the same test strip.

Turn the meter off and remove the test strip. Re-read the instructions about applying blood (starting on page 41) and repeat the test with a new test strip.

Display	Description of error
	<p><b>Error: measurement</b> Measurement error</p>
	<p><b>Error: quality control</b> The test strip failed the internal quality control check. The test strip is unusable.</p>


**Solution**

Turn the meter off and remove the test strip. Repeat the test with a new test strip. Do not touch or remove the test strip when a test is in progress.

Turn the meter off and remove the test strip. Repeat the test with a new test strip. In case the error occurs again call Roche customer service.

**Display**

**Description of error**

**Error: measurement**

Measurement error resulting from the blood sample.

## Solution

Turn the meter off and remove the test strip. Repeat the measurement using a new test strip and blood taken from a new puncture site at the tip of another finger. Do not touch or remove the test strip when a test is in progress.

**Note on Error 7 “Measurement”:** The CoaguChek XS PT test strip may be used for patients under a combination therapy of oral anticoagulants plus heparin injections. For maximum heparin concentrations which do not interfere with the test, please refer to the package insert. Under no circumstances, however, should heparinized capillary tubes be used for sample application. If capillaries are used, please use only the dedicated CoaguChek capillary tubes. Be sure to apply the blood drop to the test strip within 15 seconds of lancing the fingertip.

### WARNING


---

This error message (7) can occur in rare cases, i.e. in patients with long coagulation times ( $> 8$  INR,  $< 5\%$  Quick). If this error message appears again when you repeat the test, your results must be checked using another method. Please contact the doctor treating you immediately.

---


**Display**

**Description of error**

**Error: internal error**

Error during the internal diagnostic test which the meter performs.

**Error: internal error**

Error during the internal diagnostic test which the meter performs.

### Solution

Turn the meter off and remove the batteries. Wait at least one minute before re-inserting the batteries in the battery compartment, and then set the date and time as described starting on page 24.

Repeat the test. If the error message persists, the meter has a fault. Please call your local customer support and service centre (see page 107).

Turn the meter off and wait for at least 2 minutes. Ensure that when you perform a test the ambient temperature is stable.

Repeat the test. If the error message persists or re-appears in the future, the meter has a fault. This error indicates possible damage to the heater plate contacts. Call your local customer support and service centre (see page 107).

## Error message during infrared downloads

Display	Description of error
 A screenshot of a device's LCD display. The display shows the time '11:23' in the top left corner and the date '15-05-15' in the top right corner. In the center, there is a graphic of an infrared remote control with several curved lines radiating from its front, representing infrared light. Below the remote control, the word 'error' is displayed in a simple, sans-serif font.	<p><b>Error: infrared interface</b> Error during communication via the infrared interface caused by stray light.</p>


**Solution**

Data downloads are performed via the infrared interface and can be disturbed by bright surrounding light, especially if the light shines directly into the infrared window on the meter or at the receiving end (e.g. on a computer).

Keep away from sources of stray light when downloading results.

## Display Information and Symbols

The meter performs an automatic display check when it is turned on. In case you suspect the display or parts of it to be defective (e. g. if you have dropped the meter), you must check by means of the display test if the meter display works properly. If your display is not functioning properly, do not perform further tests as results may be misread if a segment is missing. Contact your local customer support and service centre in this case.


If the display check is too quick, you can put it “on hold” by pressing and holding down the On/Off button **ⓘ** the next time you turn the meter on. For as long as you keep the button pressed, the display remains frozen.

These symbols may appear on the display. They have the meanings shown:

Symbol	Meaning
	Test strip (without sample application area)
	Test strip (with sample application area)
	Apply sample
	User must wait until the meter has completed an action.
24 h	24-hour time format
12 h	12-hour time format
am	Time between midnight and noon (in 12-hour time format)
pm	Time between noon and midnight (in 12-hour time format)
%Q	Results are displayed as a Quick percentage value.
Sec	Results are displayed in seconds.
INR	Results are displayed in INR units.

## Display Information and Symbols

Symbol	Meaning
	Beep tone is turned on
	Beep tone is turned off
	Result is above the chosen therapeutic range (only with INR as unit)  Result is below the chosen therapeutic range (only with INR as unit)
	<ul style="list-style-type: none"> <li>■ Results in the chosen unit</li> <li>■ Error numbers</li> <li>■ Code number</li> </ul>
	Memory contains no results or no further results
	Time in the format hour : minutes.
	Date in the format day-month-year, month-day-year, or year-month-day
	The meter is in Setup Mode.
	Gives the code number of the code chip presently inserted.

Symbol	Meaning
	Result in the chosen unit is above the measuring range.
	Result in the chosen unit is below the measuring range.
	<p>Battery status:</p> <ul style="list-style-type: none"> <li>▪ When the batteries still have their full charge, all segments are lit.</li> <li>▪ Individual segments disappear one by one as the batteries become weaker.</li> <li>▪ When there is no segment remaining, you can no longer perform a test. You can, however, still access the meter's memory.</li> </ul>
<b>mem</b>	The meter is in Memory Mode.
	Automatic quality control completed successfully
<b>C</b>	Result for this blood test is flagged with a "c" due to a low haematocrit value or wet hands.
<b>error</b>	Reports an error (see Error messages)
	Room or meter temperature is outside the acceptable range.
	Communication is taking place via the infrared interface

<b>Symbol</b>	<b>Meaning</b>
	Test strip guide cover is open.

## Further Information

### Ordering

Please contact your pharmacist or specialist supplier.

---

CoaguChek XS PT Test, 24 test strips	REF 0 4625358
---	---------------

---

CoaguChek XS PT Test, 2 x 24 test strips	REF 0 4625315
---	---------------

---

CoaguChek XS Softclix Lancet, 50 lancets	REF 0 3506509
---	---------------

---

### Product limitations

Please read the information in the package insert supplied with the test strips for detailed product data and limitations.

## Product Specifications

### Operating conditions and technical data

Temperature range	+15 °C to +32 °C
Relative humidity	10 to 85% (no condensation)
Maximum altitude	4300 m
Position	Place the meter on a level, vibration-free surface or hold it so it is roughly horizontal.
Measuring range	%Q: 120 – 5 SEC: 9.6 – 96 INR: 0.8 – 8.0
Memory	300 test results with date and time
Interface	Infrared interface, LED/IRED Class 1
Battery operation	4 x 1.5 V batteries (alkali-manganese) type AAA (LR03)
Number of tests per set of batteries	Up to 2 years or up to 300 tests (depending on frequency of use)
Safety class	III
Automatic power-off	After 3 minutes
Dimensions	138 x 78 x 28 mm
Weight	127 g (without batteries)

**Sample material**

Sample type	Capillary whole blood or non-anticoagulated, venous whole blood.
Sample size	At least 8 $\mu\text{L}$
Interactions	Refer to test strip package insert

**Storage and transport conditions**

Temperature range	-25 °C to +70 °C
Relative humidity	10 to 85% (no condensation)

### **Disposal of the CoaguChek XS meter**

During measurements the meter itself may come into contact with blood. Used meters therefore carry a risk of infection. Please dispose of your used meter – after removing the batteries and cleaning the meter – according to the regulations applicable in your country. For information about correct disposal please contact your local council.

The meter falls outside the scope of the European Directive 2002/96/EC (Directive on waste electrical and electronic equipment (WEEE)).

### **For healthcare professionals**

Dispose of the meter in line with your facility's procedures for dealing with potentially biohazardous waste.

### **Disposal of used batteries**

Do not dispose of the batteries with normal domestic waste. Dispose of used batteries in an environmentally responsible manner and in accordance with applicable local regulations and directives. Contact your local council / local authority or the manufacturer of the used batteries for advice on correct disposal.

## Information service

If you have any further questions, please contact your local customer support and service centre.

---

### United Kingdom

Roche Diagnostics Ltd.  
Charles Avenue  
Burgess Hill, RH15 9RY

### CoaguChek Careline

0808 100 7666

### Australia

Roche Diagnostics Australia Pty Ltd.  
ABN 29 003 001 205  
31 Victoria Avenue  
Castle Hill, NSW, 2154  
Telephone 02-9860 2222 or  
1800 645 619

---

### Canada

Roche Diagnostics  
201 Boulevard Armand-Frappier  
Laval, Québec H7V 4A2

### Coagulation Info-Line

Telephone 1-877-426-2482

---

### New Zealand

Roche Diagnostics N.Z. Ltd  
15 Rakino Way, Box 62-089  
Mt. Wellington, Auckland  
Telephone +64 9 276 4157  
Toll Free 0800 652 634

---

## Repairs

Please note that repairs and other modifications to the meter may only be performed by persons authorized by Roche Diagnostics.

## Guarantee

The statutory guarantee provisions on rights in consumer goods sales in the country of purchase shall apply.

ACCU-CHEK, COAGUCHEK, SAFE-T-PRO and SOFTCLIX are trademarks of Roche.

© 2015 Roche Diagnostics. All rights reserved.

## Alphabetical Index

### A

Altitude (operating conditions) .....	11
Applying blood .....	52

### B

Batteries	
inserting .....	18–20
types .....	16, 18
Beep tone .....	50
setting .....	35
Blood sample (testing) .....	41–62

### C

Capillary blood sample (testing) .....	41–62
Cleaning .....	69–75
test strip guide .....	73
Coagulation process .....	53
Coagulation units .....	34, 55, 58
Code chip .....	21–23, 49
code number .....	22, 49
inserting .....	22
Collecting blood .....	44

## Alphabetical Index

<b>D</b>	
Date .....	30
Date format .....	29
Disinfecting .....	69–75
Display	
function check .....	19, 46
symbols .....	99–102
Disposal .....	16, 18, 22, 23, 56, 106
<b>E</b>	
Electromagnetic interference .....	11
Error messages .....	77–97
<b>F</b>	
Fingertip .....	51
<b>G</b>	
Guarantee .....	108
<b>H</b>	
Humidity .....	11
<b>I</b>	
Information service .....	107
Infrared interface .....	65, 96
<b>L</b>	
Lancing (your finger) .....	51

**M**

M button .....	27, 63
Measuring range .....	59
Memory .....	63–68
erasing .....	66–68
Meter	
setup .....	24–39
Meter guide .....	15

**O**

Operating conditions .....	11, 104
Ordering .....	103

**P**

Prothrombin time .....	9
------------------------	---

**Q**

Quality control .....	54
-----------------------	----

**R**

Result .....	55
selecting the coagulation unit .....	34

## Alphabetical Index

<b>S</b>	
Set button .....	26
Setup .....	24–39
beep tone .....	35
date .....	30
date format .....	29
setting up the meter .....	26
summary .....	25
therapeutic range .....	36–37
time .....	33
time format .....	32
unit .....	34
Starting up .....	17–23
Stray light .....	65, 96
Symbols	
blood drop .....	50
clr .....	67
error .....	77–97
hourglass .....	50, 55, 68
mem .....	63
on packaging .....	2
summary .....	99–102
type plate .....	2

**T**

Target range .....	36, 57
setting .....	36–37
Technical data .....	104
Telephone number (information service) .....	107
Temperature range .....	11, 42
Test	
performing .....	46–56
preparing .....	45
principle .....	10
Test strip	
inserting .....	48
symbol .....	47
Test strip guide cover .....	73
Therapeutic target range .....	36, 57
setting .....	36–37
Time .....	33
Time format .....	32

**U**

Unit	
changing temporarily .....	58, 65
results (INR, %Quick, seconds) .....	34, 55, 58

This page intentionally left blank.

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Roche Diagnostics GmbH