
HemoCue Hb 201+
Operating Manual
Manuel d'utilisation

HemoCue Hb 201+ *Operating Manual*

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
US

The HemoCue Hb 201+ Analyzer provides determination of hemoglobin quickly, easily and with lab quality results. Capillary, venous or arterial whole blood may be used. This manual provides the basic instructions for use as well as technical specifications. Additional information may be obtained from HemoCue, Inc.

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Components
Éléments

US

1. The HemoCue Hb 201⁺ Analyzer*
2. AC adapter**
3. 4 type AA batteries***
4. A vial of HemoCue Hb 201 Microcuvettes***
Individually packaged HemoCue Hb 201 Microcuvettes***
5. HemoCue Hb 201⁺ Operating Manual, HemoCue Hb 201⁺ Quick Reference Guide and Instruction CD
6. HemoCue Cleaner

The HemoCue Hb 201⁺ and its adapter are delivered in a carton. Open the carton on a stable surface and lift out the instrument and accessories. On the bottom of the analyzer there is a lid covering the battery compartment. Press the flap to remove the lid. Place 4 type AA batteries in the battery compartment, observing the indication of polarity in the battery holders. Replace the lid. Consult local environmental authorities for proper disposal of batteries.

* Do not open the cover of the analyzer.

**  Only use adapters listed under specifications.

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
Start-up
Mise en service

US

- 1a. If AC power is available, plug the supplied power adapter into the power inlet at the back of the analyzer.
- 1b. If no power is available, insert the 4 type AA batteries, 1.5 V, into the battery compartment.
- 1c. If a battery symbol appears on the display, the batteries are running low on power. The analyzer will continue to give accurate results, but the batteries should be replaced as soon as possible.
2. Pull the cuvette holder out to the loading position.
3. Press and hold the left button until the display is activated (all symbols appear on the display).
4. The display shows the version number of the program, after which it will show "⌘" and "Hb" and "🔊" (audio signal). During this time the analyzer will automatically verify the performance of the optronic unit by performing an automatic SELFTEST.
5. After 10 seconds, the display will show three flashing dashes and the HemoCue symbol. This indicates that the HemoCue Hb 201⁺ has passed the SELFTEST and is ready for use. If the SELFTEST fails, an error code will be displayed.

To turn the analyzer off, press and hold the left button until the display reads OFF and then goes blank. Disconnect the equipment from the power source by removing the AC adapter from the outlet.

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
Mesure Sang capillaire




US

1. To perform a test using capillary blood, the cuvette holder should be in its loading position. The display will show three flashing dashes and the HemoCue symbol.
2. Make sure the patient's hand is warm and relaxed. Use only the middle or ring finger for sampling. Avoid fingers with rings on.
3. Clean the finger with alcohol or a suitable disinfectant and allow to dry or wipe off with a dry, lint-free wipe.
4. Using your thumb, lightly press the finger from the top of the knuckle towards the tip. This stimulates the blood flow towards the sampling point.

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


US

5. For best blood flow and least pain, sample at the side of the fingertip, not in the center.
6. While applying light pressure towards the fingertip, puncture the finger using a lancet.
7. Wipe away the first 2 or 3 drops of blood.
8. Re-apply light pressure towards the fingertip until another drop of blood appears.

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
US

9. When the blood drop is large enough, fill the microcuvette in one continuous process. Do NOT refill!
10. Wipe off excess blood from the outside of the microcuvette with a clean, lint-free wipe, being careful not to touch the open end of the microcuvette, which could result in blood being drawn out of the microcuvette.
11. Look for air bubbles in the filled microcuvette. If present, discard the microcuvette and fill a new microcuvette from a new drop of blood. Small bubbles around the edge can be ignored.

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US

12. Place the filled microcuvette in the cuvette holder. This must be performed within ten minutes after filling the microcuvette!
13. Gently slide the cuvette holder to the measuring position.
14. During the measurement "g" and three fixed dashes will be shown on the display.
15. After 15–60 seconds, the hemoglobin value of the sample is displayed. The result will remain on the display as long as the cuvette holder is in the measuring position. When operating on battery power, the analyzer will automatically turn off after approximately 5 minutes.
16. Although the reagents are present in the microcuvette in extremely low quantities, consult local environmental authorities for proper disposal. Always handle blood specimens with care, as they might be infectious.

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

Measuring Control material, venous or arterial blood

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


④

US

The HemoCue Hb 201 system can be verified on the days of testing through the use of commercially available hematology controls. Contact HemoCue, Inc. for control information.

1. The blood or control material should be well mixed prior to performing a measurement.
- 2ab. Place a drop of blood or control material onto a hydrophobic surface, e.g., a plastic film, using a pipette or a suitable blood transfer device.
3. Fill the microcuvette in one continuous process. Do NOT refill! Wipe off excess blood from the outside of the microcuvette with a clean, lint-free wipe, being careful not to touch the open end of the microcuvette, which could result in blood being drawn out of the microcuvette.
4. Look for air bubbles in the filled microcuvette. If present, discard the microcuvette and fill a new microcuvette from a second drop of sample. Small bubbles around the edge can be ignored.

Perform the analysis as per steps 12–16 on page 15.

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
Maintenance Entretien

US

The cuvette holder should be cleaned after each day of use.

1. Check that the analyzer is turned off. The display should be blank.
2. Pull the cuvette holder out to its loading position. Carefully press the small catch positioned in the upper right corner of the cuvette holder.
3. While pressing the catch, carefully rotate the cuvette holder towards the left as far as possible. Carefully pull the cuvette holder away from the analyzer.
4. Clean the cuvette holder with alcohol or mild detergent. It is important that the cuvette holder is completely dry before being replaced.

The optronic unit should be cleaned when directed to do so in the Troubleshooting Guide or as desired.

A dirty optronic unit may cause the analyzer to display an error code. To clean the optronic unit, proceed as follows.

5. Push the HemoCue Cleaner swab into the opening of the cuvette holder. Move from side to side 5–10 times. If the swab is stained, repeat with a new swab. No further cleaning is required if the swab remains clean.
6. Wait 15 minutes before replacing the cuvette holder and using the analyzer.

The cover may be cleaned with alcohol or a mild soap solution.

As an alternative to the HemoCue Cleaner, a cotton tip swab moistened with a alcohol (without additive) or water may also be used for cleaning.

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US Trouble Shooting Guide

If you are unable to resolve the problem by following this Trouble Shooting Guide, please contact HemoCue Inc. The analyzer has no serviceable parts.

Symptom	Explanation	Action
The analyzer shows an error code.	May be a temporary fault.	Turn off the analyzer and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement. If the problem continues, see specific error code below.
E00	No stable endpoint is found within the time range. 1. The cuvette is faulty. 2. The circuit board is out of order.	1a. Check the expiration date for the microcuvettes. 1b. Take a new microcuvette and repeat the measurement. 2. The analyzer needs service. Contact HemoCue, Inc.
E01–E05	1. Dirty optronic unit or faulty electronic or optronic unit.	1a. Turn off the analyzer and clean the optronic unit as described in the maintenance section. 1b. The analyzer needs service. Contact HemoCue, Inc.
E06	1. Unstable blank value The analyzer might be cold.	1. Turn off the analyzer and allow it to reach room temperature. If the problem continues, the analyzer needs service. Contact HemoCue, Inc.
E07	1. The battery power is too low.	1a. The batteries need to be replaced. Turn off the analyzer and replace the batteries, 4 type AA. 1b. Use the power adapter.
E08	The absorbance is too high. 1. An item is blocking the light in the cuvette holder.	1a. Check that the analyzer and microcuvettes are being used according to the HemoCue Hb 201+ operating manual and instructions for use. 1b. The analyzer needs service. Contact HemoCue, Inc.
E09–E30	1. Dirty optronic unit or faulty electronic or optronic unit.	1a. Turn off the analyzer and clean the optronic unit as described in the maintenance section. 1b. The analyzer needs service. Contact HemoCue, Inc.


Symptom	Explanation	Action
HHH	<ol style="list-style-type: none"> 1. Measured value exceeds 25.6 g/dL (256 g/L, 15.9 mmol/L). 	
No characters on the display	<ol style="list-style-type: none"> 1. The analyzer is not receiving power. 2. If on battery power, the batteries need to be replaced. 3. The display is out of order. 	<ol style="list-style-type: none"> 1a. Check that the power adapter is connected to the AC power supply. 1b. Check that the power adapter is securely connected to the analyzer. 1c. Check that the cable is not damaged. 2. Turn off the analyzer and replace the batteries, 4 type AA. 3. The analyzer needs service. Contact HemoCue, Inc.
The display gives erroneous characters.	<ol style="list-style-type: none"> 1. The display is out of order. 2. The microprocessor is out of order. 	<ol style="list-style-type: none"> 1. The analyzer needs service. Contact HemoCue, Inc. 2. The analyzer needs service. Contact HemoCue, Inc.
The display shows "☐" .	<ol style="list-style-type: none"> 1. The batteries need to be replaced. 2. If on mains power, the mains adapter or the circuit board is out of order. 	<ol style="list-style-type: none"> 1. Turn off the analyzer and replace the batteries, 4 type AA. 2a. Check that the power adapter is properly connected and working. 2b. The analyzer needs service. Contact HemoCue, Inc.
The display does not switch from "2" and "Hb" to three flashing dashes and "☉" (ready for measuring).	<ol style="list-style-type: none"> 1. The magnet in the cuvette holder may be missing. 2. The magnetic sensor is out of order. 	<ol style="list-style-type: none"> 1. The analyzer needs service. Contact HemoCue, Inc. 2. The analyzer needs service. Contact HemoCue, Inc.

Symptom	Explanation	Action
<p>Measurements on control materials are out of range – either too HIGH or too LOW.</p>	<ol style="list-style-type: none"> 1. The microcuvettes are beyond their expiration date, damaged or have been improperly stored. 2. The optical eye of the microcuvette is contaminated. 3. The control has not been mixed properly and/or is not at room temperature. 4. Air bubbles in the microcuvette. 5. The optronic unit is dirty. 6. The control is not suitable for use with the HemoCue Hb 201⁺ system. 7. The calibration of the analyzer has been changed. 8. The controls are beyond their expiration dates or have been improperly stored. 	<ol style="list-style-type: none"> 1. Check the expiration date and the storage conditions of the microcuvettes. 2. Remeasure the sample with a new microcuvette. 3. Make sure that the control is mixed properly and at room temperature. 4. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 5. Clean the optronic unit as described in the maintenance section. 6. Only use controls intended for the HemoCue Hb 201⁺ systems. Contact HemoCue, Inc. for control information. 7. The analyzer needs service. Contact HemoCue, Inc. 8. Check the expiration date and the storage conditions of the control. Take a new microcuvette and repeat the measurement from a new vial/bottle of control.
<p>Measurements on patient samples are higher or lower than anticipated.</p>	<ol style="list-style-type: none"> 1. Improper sampling technique. 2. The microcuvettes are beyond their expiration date, damaged or have been improperly stored. 3. The optical eye of the microcuvette is contaminated. 4. Air bubbles in the microcuvette. 5. The optronic unit is dirty. 6. The calibration of the analyzer has changed. 	<ol style="list-style-type: none"> 1. See pages 8–17 in this manual. 2. Check the expiration date and the storage conditions of the microcuvettes. Check the entire system with a commercial control. 3. Remeasure the sample with a new microcuvette. 4. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 5. Clean the optronic unit as described in the maintenance section. 6. The analyzer needs service. Contact HemoCue, Inc.

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Specifications

General

The HemoCue Hb 201⁺ is a system used for the determination of the total amount of hemoglobin in whole blood. The system consists of a specially designed analyzer with specially designed microcuvettes containing dried reagents. The microcuvette serves as pipette, reaction vessel and as a measuring microcuvette. No dilution is required. The hemoglobin measurement takes place in the analyzer, which follows the progress of the reaction until the steady state has been reached. The system is factory calibrated against the hemoglobincyanide (HiCN) method, the international reference method for the determination of the hemoglobin concentration in blood, and needs no further calibration.

Intended Use

Quantitative determination of hemoglobin in capillary, venous or arterial blood. The HemoCue Hb 201⁺ system is used for quantitative determination of hemoglobin in whole blood using a specially designed analyzer, the HemoCue Hb 201⁺ Analyzer, and specially designed microcuvettes, the HemoCue Hb 201 Microcuvettes. HemoCue Hb 201 Microcuvettes are for *In Vitro* Diagnostic use only. The HemoCue Hb 201⁺ Analyzer is only to be used with HemoCue Hb 201 Microcuvettes.

Theory

The reaction in the microcuvette is a modified azidemethemoglobin reaction. The erythrocyte membranes are disintegrated by sodium deoxycholate, releasing the hemoglobin. Sodium nitrite converts the hemoglobin iron from the ferrous to the ferric state to form methemoglobin, which then combines with azide to form azidemethemoglobin.

Storage and handling

Use the HemoCue Hb 201 Microcuvettes prior to their expiration date. The expiration date is printed on each package.

Storage of microcuvettes kept in a vial

The microcuvettes are to be stored at room temperature (59–86 °F, 15–30 °C). Do not refrigerate. The expiration date of the HemoCue Hb Microcuvettes in a sealed vial is printed on the vial.

Once the seal is broken, the microcuvettes are stable for three months. Always keep the container properly closed.

Storage of the individually packaged cuvettes

The microcuvettes are to be stored at room temperature (59–86 °F, 15–30 °C). Do not refrigerate. The individual packaged microcuvettes are stable until the date printed on each package.

HemoCue Hb 201⁺ Analyzer

The analyzer can be stored at 32–122 °F (0–50 °C). Operating temperature is 59–86 °F (15–30 °C). Allow the analyzer to reach ambient temperature before use. The analyzer should not be operated at high (i.e. > 90 % non-condensing) humidity.

Specimen collection and preparation

Capillary, venous or arterial blood may be used. Appropriate anticoagulants (e.g. EDTA or heparin) may be used, preferably in solid form to avoid dilutional effects. Mix all specimen tubes thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 8–10 times by hand. Hemoglobin remains unchanged for days, provided that the blood does not become infected. If the specimen has been stored in a refrigerator, it will be viscous and the blood should be allowed to warm up to room temperature before mixing².

In case of peripheral circulatory failure, hemoglobin measurements from capillary sampling may be misleading. Confirmation of an unacceptable or unexpected laboratory result is considered as a common, well established routine in the laboratory. Extra cellular fluid present in capillary samples due to individual physiology or the affects of unacceptable capillary sampling technique can cause

false hemoglobin results.

One of several pre analytical factors associated with technique revolves around the size and action mode of the lancet and how it is used. This and other factors associated with capillary sampling technique can have an impact on the puncture, blood flow and ultimately on the results obtained. Additional factors unrelated to sampling technique, such as an individual's capillary circulation status, must also be recognized as potential causes of false results. Therefore when testing for hemoglobin using capillary blood as the sample, confirmation of an unacceptable or unexpected result is needed to rule out potential pre analytical factors as the cause. This is particularly helpful when a cut off value is used, i.e. when screening for nutritional status based on hemoglobin. Since the individual performing the test may not always be aware of all the potential pre analytical factors or even those that are directly or indirectly associated with sampling techniques, a second finger stick should be performed when wanting to rule out pre analytical factors as a cause of an unacceptable or unexpected result. Performance of the second finger stick by a different individual may aid in lessening the affect of the factors that could have been technique related.

Always handle blood specimens with care, as they might be infectious. Consult local environmental authorities for proper disposal. Always wear protective gloves when handling blood specimens. The microcuvette is for single use only.

Materials required

- HemoCue Hb 201⁺ Analyzer
- HemoCue Hb 201 Microcuvettes
- Lancet (for capillary samples)
- Pipette or other transfer device (for venous, arterial or control material samples)
- Lint-free wipe (non-fraying)

Hydrophobic plastic or glass slide (for venous, arterial or control material samples)

Quality Control

The HemoCue Hb 201⁺ Analyzer has an internal quality control, the "selftest". Every time the analyzer is turned on, it will automatically verify the performance of the optronic unit of the analyzer. This test is performed every second hour if the analyzer remains switched on.

Upon passing the SELFTEST, the display will show the HemoCue symbol and three flashing dashes, indicating that the analyzer is ready to perform a measurement. An error code will be displayed if the SELFTEST fails. If use of liquid control material is required by local or other regulations, contact your local distributor for control information.

Expected values²

Adult Males 13.0–17.0 g/dL (130–170 g/L, 8.1–10.5 mmol/L)

Adult Females 12.0–15.0 g/dL (120–150 g/L, 7.4–9.3 mmol/L)

Infants, after neonatal period 11.0–14.0 g/dL (110–140 g/L, 6.8–8.7 mmol/L)

Children, two years to teenage: gradual increase to adult normals. Due to a wide range of conditions (dietary, geographical, etc) which affect normal values, it is recommended that each laboratory establish its own normal range.

Measuring range

0–25.6 g/dL (0–256 g/L, 0–15.9 mmol/L).

Results above 25.6 g/dL (256 g/L, 15.9 mmol/L) will be displayed as HHH. Values above 23.5 g/dL (235 g/L, 14.6 mmol/L) must be confirmed using a suitable laboratory method.

Limitations

HemoCue Hb 201 Microcuvettes are for *in vitro* diagnostic use only. The HemoCue Hb 201⁺ Analyzer is only to be used together

with HemoCue Hb 201 Microcuvettes. For further limitations of the procedure, see the HemoCue Hb 201 Microcuvettes package insert.

Technical Specifications

Dimensions: 3.35 × 6.30 × 1.69 inches (853 × 160 × 343 mm).

Weight: 0.77 pounds (350 g) (with batteries installed).

4 batteries type AA.

Power adaptor: CE marked.

Only use adapters recommended by HemoCue or equivalent, as listed on page 34.

Pollution degree: 2.

Overvoltage category: II.

Atmospheric pressure: 700 hPa to 1060 hPa.

Equipment not suitable for use in the presence of flammable mixtures.

The instrument is tested according to IEC 61010-1, First edition, 1990 with Amendment 1, 1992 and Amendment 2, 1995 (EN 61010-1: 1993 with A2: 1995), EN 60601-1-2 and complies with the IVD Medical Device Directive 98/79/EC. Signal input and signal output

Accessory equipment to the analog and digital interfaces must be certified according to the respective UL standards (i.e. UL 60950-1 for data processing equipment and IEC 6060-1 for medical equipment). Furthermore all configurations shall comply with the system standard IEC 6060-1. Connections of additional equipment to the signal input or signal output connections is deemed to be “configuring a medical system”, and therefore assumes responsibility for the systems compliance within the requirements of IEC 60601-1-1. If in doubt, contact HemoCue AB.

Warranty

The analyzer carries a 24-month warranty from the day of receipt. After the warranty period, service/repair is carried out at fixed prices.

Spare parts – Accessories

The following spare parts and accessories are available:

Power adaptor

Battery Lid

HemoCue Cleaners

Hard Carrying Case

Cuvette holder

HemoCue Lancets

HemoCue Hb 201⁺ Control Cuvette

Recommended separation distance between Portable and mobile RF communications equipment and HemoCue Hb 201+ Analyzer

The HemoCue systems are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of HemoCue systems can help prevent electromagnetic interference by maintaining a minimum distance between portable and RF communications equipment (transmitters) and HemoCue systems as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at maximum output power not listed above, the recommended separation distances (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where (P) is the maximum output power rating of the transmitter in watts (W) according to transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – Electromagnetic immunity

The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or user of the HemoCue systems should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	< 5 %U (> 95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles < 5 %U (> 95 % dip in U) for 5 seconds For explanation of U see NOTE 1	< 5 %U (> 95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles < 5 %U (> 95 % dip in U) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the HemoCue systems requires continued operation during power mains interruptions, it is recommended that the HemoCue systems be powered from an uninterruptible power supply or a battery.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p data-bbox="302 254 540 306">Conducted RF IEC 61000-4-6</p> <p data-bbox="302 336 540 389">Radiated RF IEC 61000-4-3</p>	<p data-bbox="540 254 777 306">3 Vrms 150 kHz to 80 MHz</p> <p data-bbox="540 336 777 389">3 V/m 80 MHz to 2.5 GHz</p> <p data-bbox="540 418 777 442">See NOTE 2 and NOTE 3</p>	<p data-bbox="777 254 950 283">3 Vrms</p> <p data-bbox="777 336 950 365">3 V/m</p>	<p data-bbox="950 254 1347 412">Portable and mobile RF communications equipment should be used no closer to any part of the HemoCue systems, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p data-bbox="950 418 1347 442">Recommended separation distance:</p> <p data-bbox="950 471 1347 495">$d=1.2\sqrt{P}$</p> <p data-bbox="950 506 1347 530">$d=1.2\sqrt{P}$ 80 MHz to 800 MHz</p> <p data-bbox="950 542 1347 565">$d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p data-bbox="950 595 1347 700">Where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and (d) is the recommended separation distance in meters (m).</p> <p data-bbox="950 730 1347 836">Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b).</p> <p data-bbox="950 865 1347 912">Interference may occur in the vicinity of equipment marked with the following symbol.</p> <div data-bbox="1117 941 1182 1006" style="text-align: center;">  </div>

NOTE 1 U is the a.c. mains voltage prior to application of the test level.

NOTE 2 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 3 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HemoCue systems are used exceeds the applicable RF compliance level above, the HemoCue systems should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the systems.
- b) Over the frequency range 150 kHz to 80 MHz, field strength should be less than 3 V/m.

Technical specifications (EMC-RF)**Use only cables with the following specification:**

USB shielded maximum 2 m

Serial shielded maximum 1.5 m

Guidance and manufacturer's declaration – electromagnetic emissions		
The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or the user of the HemoCue systems should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions	Group 1	The HemoCue systems uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment
RF emissions	Class B	The HemoCue systems are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Symbols used

Attention, see instructions for use

CE mark

Refers to Type BF Applied Part, Per EN60601-1

Class II equipment

The Medical Equipment with respect to electrical shock, fire and mechanical hazards only in accordance with UL-2601-1 and related collateral IEC 60601-1 standards for which the products has been found to comply by UL.

Only valid within the European Community. Indicates separate collection for waste of electrical and electronic equipment.

Temperature limitation

Efficiency Level

References

1. Makarem, A. In clinical Chemistry: Principles and Technics, 2nd ed., Henry, R.J., Cannon, D.C., and Winkelman, J.W., Harper and Row, Hagerstown, M.D., 1974, pp. 1125–1147.
2. Dacie and Lewis, Practical Haematology Ninth edition.
3. Jacques Wallach, M.D., Interpretation of Diagnostic Tests.

Manufacturer

HemoCue AB
 Box 1204
 SE - 262 23 Ängelholm, Sweden
 Phone: +46 77 570 02 10
 Fax: +46 77 570 02 12
 E-mail: info@hemocue.se
 www.hemocue.com

HemoCue Distributor USA:

HemoCue, Inc.
 40 Empire Drive
 Lake Forest, CA 92630
 Phone (General) 800-881-1611
 Orders 800-323-1674
 Technical support 800-426-7256
 Fax (Cust. Service) 800-333-7043

Recommended Adapter or equivalent:

Country	Type	Input
USA/Canada	Mod nr FE3515060D035	120 V~/60 Hz/6.6 VA
World wide	HCA01	100 V~/240 V~/50–60 Hz, < 500 mA

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HemoCue AB, Box 1204, SE-262 23 Ängelholm, Sweden. Phone +46 77 570 02 10. Fax +46 77 570 02 12.