

HOW TO TRIM YOUR EARMOULD TUBING



•Step 1: Your ear-mould will arrive with long tubing, like this

•Step 2: Put the ear-mould in your ear, ensuring you insert it correctly.



•Step 3: Put the hearing aid behind the ear and mark the point where the hearing aid meets the tubing (see arrow)



•Step 4: Trim the tubing about quarter of an inch/3mm up from the marked point and then connect the tubing to the hearing aid elbow (hook)

Audiology Department (C41)
Royal United Hospital
Combe Park
Bath
BA1 3NG

01225 824035

ruh-tr.audiology@nhs.net